
SENCOR 

*Wrist blood
pressure monitor*

SBP 100

Instruction manual EN

Blood pressure Monitor

Instruction manual

SBP 100 is a blood pressure monitor that maximizes a reading accuracy of blood pressure through the artery in your wrist echoed from the heart. It detects arrhythmia symptoms and displays the symbol „“ on LCD. If arrhythmia symptoms symbol is displayed, you are recommended to consult with a doctor. Please do not interpret Blood pressure measurements including arrhythmia symptoms by yourselves. The reading should be interpreted by a physician or a trained health professional who is familiar with your medical history. By using the device regularly, and recording the results for your physician to interpret, you keep your physician informed of the continuing trends in your blood pressure.

The measuring accuracy of SBP 100 was tested at the time of manufacture and was proven by clinical research in accordance with AAMI/ANSI-SP10.

Blood pressure Monitor


Instruction manual

EN

What you should know about blood pressure

What is Blood Pressure ?

Blood pressure is the force of blood flowing against the walls of the arteries. Arterial blood pressure varies during the course of the cardiac cycle. The highest pressure in the cycle is called the **SYSTOLIC BLOOD PRESSURE** ; the lowest the **DIASTOLIC BLOOD PRESSURE**. Both pressure readings, the **SYSTOLIC** and **DIASTOLIC**, are necessary to enable a physician to evaluate the status of a patient's blood pressure. Many factors such as physical activity, anxiety, or the time of day, can influence your blood pressure. Blood pressure constantly changes throughout the day. It rises sharply in the early morning and declines during the late morning. Blood pressure rises again in the afternoon and finally drops to a low level at night. Also, it may vary in a short period of time. Therefore, reading from successive measurements can be fluctuate.

Blood pressure Monitor

Instruction manual

Why is it important to measure your Blood Pressure at home?

Blood pressure measured in a doctor's office only provides a momentary value. Repeated measurements at home better reflect one's actual blood pressure values under everyday conditions. Moreover, many people have a different pressure when they measure at home, because they tend to be more relaxed than when in the doctor's office. Regular blood pressure measurements taken at home can provide your doctor with valuable information on your normal blood pressure values under actual "everyday" conditions.

The World Health Organization(WHO) Blood Pressure Classifications

WHO has set up the following standard blood pressure values when measured at resting pulse.

Blood	Normal values	Mild hypertension	Severe hypertension
Pressure(mmHg)			
SYS =Systole = upper value	up to 140	140 – 180	over 180
DIA =Diastole= lower value	up to 90	90 – 110	over 110

What are arrhythmias?

The heart is a pump that functions by pushing the blood through its four chambers. The blood is „pushed“ through in a controlled sequence of muscular contractions. The sequence is controlled by bundles of cells which control the electrical activity of the heart. When the sequence is disturbed, heart arrhythmias occur.

Blood pressure Monitor

Instruction manual

EN

Arrhythmias are abnormal rhythms of the heart. They cause the heart to pump blood less effectively. Most cardiac arrhythmias are temporary and benign. Most temporary and benign arrhythmias are those where your heart skips a beat or has an extra beat. These occasional changes can be brought on by strong emotions or exercise. Nonetheless, some arrhythmias may be life-threatening and require treatment.

Symptoms?

Symptoms of Arrhythmia overall: Palpitations or rapid thumping in your chest, feeling tired or light-headed, losing consciousness, shortness of breath and chest pain.

Symptoms of bradycardia: feeling tired, short of breath, dizzy, or faint.

Symptoms of tachycardia: heartbeat might feel like a strong pulse in your neck, or a fluttering, racing beat in your chest, feelings of discomfort, weak, shortness breath, faint, sweaty, and dizzy.

Can arrhythmia be treated?

Treatment of arrhythmia depend on the type of arrhythmia, the patients age, physical condition and age. Methods are available for prevention of arrhythmia. These methods include relaxation techniques to reduce stress, limit intake of caffeine, nicotine, alcohol and stimulant drugs. Many arrhythmias require no treatment, they are naturally controlled by the body's immune system. However if it is necessary that arrhythmias must be controlled, they can be controlled by drugs, cardioversion, automatic implantable defibrillators or an artificial pacemaker. Arrhythmias are very serious. They cause nearly 250,000 deaths in the United States each year. Arrhythmias that start in the lower chambers of the heart (the ventricles) are more serious than those that start in the upper chambers (the atria).

Blood pressure Monitor

Instruction manual

Arrhythmia statistics

- During a 24-hour period about 20% of healthy adults are likely to have multiple types of premature ventricular beats
- As many as 2 million Americans are living with arterial fibrillation, a type of arrhythmia

Caution

- To ensure accurate measuring results, carefully read these use instructions.
- This product is intended for household use only. Keep product and batteries away from children.
- People suffering from cardiac arrhythmia, vascular constriction, arteriosclerosis in extremities, diabetes or users of cardiac pacemakers should consult their doctor before measuring their blood pressure themselves, since deviations in blood pressure values may occur in such cases.
- If you are under medical treatment or taking any medication, please consult your doctor first.
- This device is designed to make automatically rapid exhaust if cuff pressure is 320mmHg or more. If the cuff doesn't make automatically rapid exhaust when cuff pressure is 320mmHg or more, please detach the cuff from the wrist.
- The use of this blood pressure monitor is not intended as a substitute for consultation with your doctor.

Blood pressure Monitor

Instruction manual

EN

Contents

1. Key rules for accurate blood pressure measurement
2. Operation guide
3. Parts identification
4. Inserting batteries
5. Attaching pressure cuff
6. Correct measurement position
7. How to measure
8. What to do if
9. Specifications

1. Key rules for accurate blood pressure measurement

1. Always take readings at the same time of day, ideally the morning, under the same conditions.
2. Do not measure within 30 minutes after smoking, coffee or tea consumption.
3. Take off wrist watch and jewelry before fitting the monitor on the measuring arm.
4. While taking a measurement, sit down, relax, keep still, especially the hand of the measuring arm, do not move or speak.
5. Hold the monitor at the heart level.
6. Wait approximately 3 minutes under take-off of wrist cuff before repeating a measurements.

Blood pressure Monitor


Instruction manual

2. Operation guide

See the steps below for your convenience.

Inserting batteries	See page 9
How to measure	See page 11
Attaching pressure cuff	See page 9
Memory recall function	See page 13

3. Parts identification

A: LCD Display
C: Time/date button
E: Battery compartment

B: Memory button
D: Start/stop button
F: Cuff

Blood pressure Monitor

Instruction manual

EN

4. Inserting batteries

1. If the Low Battery Indicator appears on the display, replace both batteries.

(Low battery indicator)

2. Inserting batteries

- Use alkaline batteries, type LR03(AAA) as supplied with this product.
- Remove the battery compartment cover and insert two batteries with correct polarity (see symbol in battery compartment).
- New batteries will provide approximately 200 measurements.

 Only discard empty batteries. They should not be disposed of in the house hold waste, but at appropriate collection sites or at your retailer.

5. Attaching pressure cuff

Fasten the monitor on the palm side of your wrist, roughly 2 cm below the palm line as shown. Make sure that the cuff fits tightly and wraps correctly & formally to get accurate reading.

Blood pressure Monitor

Instruction manual

**Note**

Please use only a cuff supplied by SENCOR.
Using other cuffs may impede measurement.

6. Correct measurement position

Measurement values vary according to your posture. Measure your blood pressure in the sitting position.

- Place your elbow on a table or other object.
- Use the storage case to position the arm so that the pressure cuff is at the same height as the heart.
Make sure that the cuff is not stuck on the storage case for smooth inflation.
- Open your hand and fingers slightly, and relax.
- Blood pressure will vary in relation to the level of the cuff from the heart. A 10cm difference will result in an 8mmHg measurement error.
- Relax your hand with the palm facing up snugly.

**Note**

Your heart is located slightly below your armpit.
For successive readings, wait 3 to 5 minutes between measurements.

- Waiting allows the engorged blood vessels to return to normal before each new measurement is taken.

Blood pressure Monitor

Instruction manual

EN

7. How to measure

Please read "Key rules for accurate blood pressure measurement" and "Correct measurement position" and "Attaching pressure cuff" before taking measurement.

Turn power on

- When Power is turned on, all display symbols appear for approximately 1 second.
- The display symbols disappear and the ready symbol(---) starts to flash. If the air exits in the cuff, the air is released automatically while the "--" symbol displays.

Automatic inflation

- The pressure cuff inflates after the "0" is displayed on the window.
- When the power button is subsequently pressed during pressurization, the unit is turned off and the cuff deflates.

When inflation is completed, the "♥" symbol will appear indicating that measurement is in progress.

- Continuous rapid step air deflation
- The pressure deflation starts by opening the solenoid valve for a short time. After pulse is detected, pressure deflation is synchronized with patient's pulse cycle and deflation is carried out.

Blood pressure Monitor


Instruction manual

1. When a pulse is detected, the pressure value moves to the lower part of the display. At this time the "♥" symbol flashes at the same rate as pulse rate.

(Pressure value)

2. At the end of the measurement, the reading is displayed. When arrhythmia symptoms is detected, at the end of the measurement, "⚡" symbol is displayed at the upper left part of the display together with the reading. If arrhythmia symptoms symbol "⚡" is displayed, you are recommended to consult with a doctor.

3. "Err" symbol appears when the accurate blood pressure could not be obtained.

Turn power off

Press the start/stop button to turn off. Otherwise the product will turn off automatically after 2 minutes.

Blood pressure Monitor

Instruction manual

EN

Setting time and date

- 1) Push the time/date button "⌚". The lower display line shows "Time" and flashing digits in the hour position. To adjust the hours, repeatedly push the "memory" button.
- 2) Push button "⌚" again: The minutes digits will flash. Again, repeatedly push the "memory" button to adjust the minutes.
- 3) Push button "⌚" again to adjust the date: The lower display line shows "Date" and flashing digits in the month position. Adjust the month by repeatedly pushing the "memory" button.
- 4) Push button "⌚" again: The day digits will flash. Again, repeatedly push the "memory" button to adjust the day.
- 5) To exit the time and date setting mode, push button "start/stop".
- 6) Otherwise the device automatically turns off after 2 minutes and activates the settings made so far.

24 hours systems with military hour

(Hour Setting)

Pressing time
button

(Minute Setting)

Pressing time
button

(Month Setting)

Pressing time
button

(Date Setting)

Memory function

After each measurement, the blood pressure monitor automatically stores the reading in memory. Please note that stored data will be erased if battery is changed.

Blood pressure Monitor

Instruction manual

Recalling stored readings

- 1) To recall stored readings, please turn off the device.
- 2) Then press the memory button. The latest reading will be displayed.
- 3) Repeatedly pressing the memory button will call up one by one all stored readings.
- 4) To exit the memory mode, turn the device off by pressing the start/stop button.
- 5) The readings of the last 60 measurements with time and date are stored in memory.

(Memory Number "18")

SBP 100: The reading of the last 60 measurements with time and date are stored in memory.
The lower display line successively shows date, time and pulse of the respective measurement.

Maintenance and Storage

- Clean the product with a soft, slightly damp cloth.
- Do not wash the main unit and disassemble the unit.
- Do not drop the unit.
- Protect it from sudden jars or shocks.
- Stains on the cuff can be carefully removed with a damp cloth.
- Do not wash or dry clean the cuff.
- Do not storage the product in a place exposed to direct sunlight or high humidity.

Calibration

This product has been calibrated at the time of manufacture. If used according to the use instructions, periodic recalibration is not required.

Blood pressure Monitor

Instruction manual

EN

8. What to do if

Problem	Reason	Solution
	Batteries are flat	Insert new batteries
Err	Measurement/ Handling error	Repeat measurement as described under "how to measurement". If problem Persists, contact a Service Centre.
Power off	Inflation above 320mmHg or after 2 minutes under non-use	Repeat measurement
---	Remained air in the cuff.	Please wait for automatically releasing air
Readings from repeated measurements differ considerably	Blood pressure is a fluctuating value. For healthy adults, deviations of 10 to 20 mmHg are possible.	None
	Cuff is not properly fitted	Ensure that the monitor is fitted as follows. - on the palm side of your wrist - approx. 2cm away from the palm line.
	Readings were not Consistently take at the heart level.	For each measurement, hold the monitor at heart level.
	Taking, coughing, laughing, moving etc. when measuring will influence the reading.	While measuring, relax, keep still, do not move or speak.

9. Specifications

Model:	SBP 100
Operation principle:	Oscillometric method, fuzzy logic
Display:	Liquid Crystal Display. 13 digits, SYS/DIA/PULSE/TIME&DATE
Measuring range:	0 - 320mmHg(cuff pressure), 40 - 199beat / min (pulse rate)
Blood pressure measurement:	20 mmHg(minimum diastolic value), 300 mmHg(maximum systolic value)
Laboratory accuracy :	± 3 mmHg(cuff pressure), $\pm 5\%$ of reading(pulse rate)
Clinical accuracy:	According to AAMI-SP10 with auscultatory reference : < 5 mmHg systematic offset, < 8 mmHg standard deviation
Inflation:	Automatic
Deflation:	Continuous valve(continuous rapid step air deflation)
Power supply:	1 . 5 V $\times$ 2 batteries, type LR 03(AAA)
Operating temperature:	+10 °C to + 40 °C
Storage temperature:	- 20 °C to + 60 °C
Humidity:	Up to 85% RH
Cuff:	Suitable for wrist circumference of 13 to 21 cm
Clinical validation:	According to AAMI-SP10 with auscultatory reference (Please note that a validation against intra-arterial measurements can lead to different results for devices which are calibrated against auscultatory reference)

If device is not used within specified temperature and humidity ranges the technical accuracy can not be ensured.

Type BF equipment

Blood pressure Monitor

Instruction manual

EN

Subject to change without notice.

This device conforms to the following standards:

DIN EN 60601 – 1: 3/96 "Medical Electrical equipment" – Part 1: General requirements for safety.

DIN EN 1060 – 1: 12/95. AMD 1 09/02 " Non-invasive sphygmomanometers" – Part 1: General requirements.

DIN EN 1060 – 3: 9/97 " Non-invasive sphygmomanometers" – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems.


This product confirms to the provisions of the EC directive 93 / 42 /
EEC (medical Device Directive).

Disposal of Used Electrical & Electronic Equipment

The meaning of the symbol on the product, its accessory or packaging indicates that this product shall not be treated as household waste. Please, dispose of this equipment at your applicable collection point for the recycling of electrical & electronic equipments waste.


In the European Union and Other European countries which there are separate collection systems for used electrical and electronic product. By ensuring the correct disposal of this product, you will help prevent potentially hazardous to the environment and to human health, which could otherwise be caused by unsuitable waste handling of this product. The recycling of materials will help conserve natural resources. Please do not therefore dispose of your old electrical and electronic equipment with your household waste. For more detailed information about recycling of this product, please contact your local city office, your household waste disposal service or the shop where you purchased the product.

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